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**DRAMINSKI 4QMAST**  
Mastitis Detector for 4 Quarters

**Manufacturer:**  
**DRAMINSKI**  
**Ul. Owocowa 17**  
**10-860 Olsztyn**  
**Poland**

**Phone: +48 89 527 11 30**  
**Fax: +48 89 527 84 44**  
**e-mail: [draminski@draminski.com](mailto:draminski@draminski.com)**

[www.draminski.com](http://www.draminski.com)

**CONTENTS**

<b>1.</b>	<b>INTRODUCTION.....</b>	<b>03</b>
<b>2.</b>	<b>DESCRIPTION.....</b>	<b>04</b>
<b>3.</b>	<b>OPERATION.....</b>	<b>05</b>
<b>4.</b>	<b>TAKING A MEASUREMENTS.....</b>	<b>07</b>
<b>5.</b>	<b>INTERPRETATION OF RESULTS.....</b>	<b>09</b>
<b>6.</b>	<b>FINAL COMMENTS.....</b>	<b>12</b>
<b>7.</b>	<b>CARE &amp; CLEANING.....</b>	<b>14</b>
<b>8.</b>	<b>BATTERY REPLACEMENT.....</b>	<b>15</b>
<b>9.</b>	<b>TECHNICAL DATA.....</b>	<b>16</b>

## INTRODUCTION

**Mastitis** - inflammation of the udder - is one of the most important and costly problems in dairy farming.

Subclinical states of mastitis are particularly dangerous. In this state, milk looks normal and the udder or quarter affected appears healthy. Subclinical mastitis can rapidly become clinical mastitis or can continue at subclinical levels, affecting milk production, milk quality, and spreading infection to other cows in the herd.

Research carried out over many years has established that the development of subclinical mastitis is accompanied by a rise in the level of salt in the milk, which immediately lowers its electrical resistance. Since the discovery of this relationship, electronic methods of detecting subclinical mastitis have become widely accepted and practical.

The relationship of milk electrical resistance and health of the udder resulted in the development of the DRAMIŃSKI MASTITIS DETECTOR by Mr. J. Damiński in 1989 with the co-operation of Polish agricultural universities.

Since then, the instrument has been widely sold in Europe, the Americas, Australia and New Zealand. Its simplicity of use, robustness, speed, low power usage and overall functionality combined with a wide-scale digital readout for interpreting severity of infections, has resulted in the success of the DRAMIŃSKI Mastitis Detector over competitor devices.

## DESCRIPTION

The DRAMIŃSKI ELECTRONIC MASTITIS DETECTOR consists of a measurement cup with electrodes, electronic unit with LCD panel, handle containing a standard 9 volt battery and a switch.

**Application of a special display showing results taken from all 4 quarters makes work more effective. The detector is easy to use and its handling comfortable. Taking measurements does not require any special knowledge nor help of any other person.**


At the bottom of a measurement cup 2 metal electrodes for measuring electrical conductivity are placed.

The device housing is drop-proof what makes its use more safe and cleaning easier.

The instrument is made of polypropylene, resistant to atmospheric conditions, most chemicals, shock, breakage, and not easily wet by milk. It is sealed and waterproof, for easy washing and cleaning.

## OPERATION

The display indicates the proper work of the device. The first check-out of the detector should be conducted with no milk in the cup, in the 'empty air', when electrodes are spaced. After pressing the button four column of dashes are displayed indicating the device readiness to work. Now the consecutive milk quarters can be measured. Blinking dashes indicate the column where the results from specific milk quarter will be shown. To take a measurement the milk is to be poured in and the switch pressed.

After the switch is pressed, the flashing dashes disappear, and result will come up in specific numbers. Next, the flashing dashes will appear in the next column. Similarly all other measurements from the next quarters will be done.

Additionally, the display indicates the exhaust state of the battery. A flashing battery symbol signals that the battery is to be replaced.

10 units is the lowest result possible whereas the measurement range is 10-990 units. After the level of 990 units has been exceeded dashes are shown on the display (the same symbol as in the case when the device tried in the 'empty air').

The measurement range of Mastitis Detector exceeds the maximum conductivity value of tested milk considerably. It means that actually it will never be exceeded.

**NOTE ! You must not use Draminski Mastitis Detector for the bulk milk.**

A minimum of 15 ml of milk (indicated by the line inside the cup) is required, and the sample must be first foremilk for accurate readings. Later milk samples will give slightly different readings.

**Before starting your measurement the following should be done:**

1. Check the work of the detector in the 'open air' ensuring that the battery is fully charged.
2. Ensure that the electrodes in the cup are clean. Contamination of electrodes with skin-oil from handling or examination while on display in stores and/or dried milk-fat can result in incorrect (high) readings. If required, wipe the electrodes with methylated spirits on a clean cloth or tissue, or similarly with kitchen detergent and rinse off.
3. Prepare a bucket full of warm water for rinsing the cup after each measurement. Additionally, another empty bucket will be useful for pouring out the milk from the cup after the measurement. To comply with the hygienic regulations you must not pour the milk out on the cowshed bedding.

### **Easy test to check whether the electrodes are clean:**

- a) prepare a solution – pour a pinch of salt into a glass of water (of room temperature) dissolving it completely.
- b) fill out the cup with the solution and read the result (it depends on the solution concentration).
- c) Pour out the solution into the glass and clean thoroughly the electrodes with a tampon moisten with the liquid for washing the dishes. Rinse the cup under the tap shaking the water remaining off.
- d) pour the solution into the cup again and read the result. If there is a remarkable drop in the result in comparison to the previous measurement it means the electrodes were coated with fat.
- e) the best is to repeat cleaning procedures until it is sure that the results do not change any more. The latter means the electrodes are now clean. Note that even touching the electrode with a finger may contaminate the electrodes.

**It is highly recommended to clean the electrodes once a week be he means of washing powder . It will remove all remaining from the cup surface.**

## TAKING A MEASUREMENTS

**ALLWAYS THE FIRST SQUIRTS OF MILK SHOULD BE MEASURED that are squirted directly to the measurement cup. Each quarter should be examined separately.**

1. Press the ON/OFF switch and place the cup under the teat A. Squirt the milk directly into the cup to fill it up (minimum 1 cm from the cup edge)
2. After about 1 second press the switch to turn the detector on – the result will be displayed. For example:

<b>370</b>	<b>---</b>
<b>---</b>	<b>---</b>

3. Discard milk into bucket, than press turn on the device, the dashes will start flashing again for the next quarter, proceed as before and for each the remaining quarters.

<b>370</b>	<b>((---))</b>
<b>---</b>	<b>---</b>

4. After you have sampled four quarters, they will all be shown in the display window.

<b>370</b>	<b>380</b>
<b>320</b>	<b>380</b>

5. Press the button once again to see which quarter may be under the threat of subclinical mastitis. Now the differences between particular quarters are displayed:

<b>10</b>	<b>0</b>
<b>60</b>	<b>0</b>

6. After the results have been viewed, the detector should be thoroughly rinsed out in order to avoid contamination of the next set of readings.
7. It is possible to examine the cow at once. You only have to press the button for a short time to turn on the detector. All previous read-outs will disappear. The detector is now ready for next measurements.

<b>((---))</b>	<b>---</b>
<b>---</b>	<b>---</b>

8. The detector turns off if the button is kept pressed longer than 3 seconds. When finished press the ON/OFF switch. If no button is pressed during 30 second time the instrument turns off automatically to conserve battery power.

## INTERPRETATION OF RESULTS

### THE RESULTS SHOULD BE EVALUATED ACCORDING TO THE FOLLOWING CRITERIA:

1. What is the numeric value of examined quarter and jaką it a typical value for a particular cow (each cow should be evaluated taking into account its age)?
2. How big are the differences between the quarters of a particular cow?

#### Ref. 1:

##### **Readings below 250 units:**

This is a clear indication of a rapid increase in the severity of infection as subclinical mastitis progresses to clinical states or the high risk of passing to subclinical states of mastitis.

##### **Readings above 300 units:**

The milk sample is of high quality and is healthy. Usually the read-outs are placed in the range 330-360 units. The incidence of subclinical mastitis is very low. Among young milk cows (1-2 lactations) the most common readings approach the vicinity of 370-400 units whereas the old cows will show readings at lower level of 300-320.

##### **Readings between 300 and 250 units:**

A progressively increasing incidence of subclinical infection as readings decrease. Due to physiological differences it is extremely difficult to define border between healthy and sick quarter. The readings at the level of 250-300 units can be taken as normal readings and quarters as healthy, especially when results in a particular cow does not show higher values. However, if the sudden drop is noticed down to 250-300 units whereas all the previous examinations gave much higher results, i.e. over 300 units, it means that this particular cow may be at risk of having subclinical masittis.

#### Ref. 2:

The difference higher than 40-50 units between the highest and the lowest result in milk quarters indicates the beginning of subclinical mastitis. It is necessary to examine the cow before each squirting systematically to see whether the non-healthy state is progressing (i.e. whether the readings drop). Additionally it is necessary to follow the hygienic regulations and undertake special care before and after squirting.

#### EXAMPLES OF RESULT INTERPRETATION:

##### young cow, 2nd lactation

all quarters gave readings over 300 units, but the difference between quarter C and A suggest that that the former may be at risk of subclinical mastitis. It is necessary to continue examining the cow.

**quarter A = 370**

**quarter B = 380**

**quarter C = 310** (the difference is 70 between maximum (qrt. B)

**quarter D = 370**

##### middle-aged cow, 5th lactation

the results suggest that one of the quarters (D) may be endangered by the subclinical mastitis. Preventive actions should be undertaken promptly and the cow put under vigilance.

**quarter A = 340**

**quarter B = 350**

**quarter C = 350**

**quarter D = 260** (result below 300, the difference 90 units to the maximum (qrt. B and C)

##### old cow, 9th lactation

one of the quarters (B) the reading was below 300 units, but it is common among the cows of that age.

**quarter A = 310**

**quarter B = 290** (result below 300 but the maximum difference only 20 units)

**quarter C = 300**

**quarter D = 300**

**The following are the factors that should be taken into account and which may influence the readings in the case of healthy cows:**

1. animal age
2. breed (the above given interpretation is valid for b/w + **Holstein-Frise** cows; for Jersey breed the readings will be higher).
3. individual characteristic of particular cow, for example fat content in milk.
4. nutrition – food content or diet change may influence the readings.
5. physiological state (for example during the ovulation the result must not be interpreted).
6. metabolic disorders, (alkalosis or acidosis can be detected by noticing drop or raise of readings among all cows in the herd)

Note! Only systematic examinations of cows enable using all potential of the Draminski Mastitis Detector.

**You must not examine the quarters with visible states of mastitis and severe changes in the milk.**

**Any proper interpretation of results will not be possible.**

## FINAL COMMENTS

1. **The point or reading below which treatment is required can not be precisely defined** and is also determined by farmer (or veterinary) choice. Quarters showing readings of 280 units or less should be regularly monitored as any rapid reduction in reading is an indication of increasing infection or of mechanical injury occurring.
2. **Each animal ought to be treated individually** therefore slight variations from typical results presented in this instruction manual do not necessarily indicate illness.  
To accurately evaluate udder health and detect subclinical mastitis simply observe the readings and variations of results typical for the given cow.  
For dairy cows giving average readings ranging between 330-340, only a decrease below 300 units is an indication of a mastitis threat. For animals with typical results in the range 370 - 380, a reading of 330 units is a similar warning.
3. **Among young milk cows (1-2 lactations)** the most common readings approach the vicinity of 400 units (ex. 2).
4. **Older cows (e.g. 9 years or older)** physiologically have increased levels of *salt and somatic cell count* in their milk and can typically result in readings of around 290 units even though they are healthy. It is important with such older cows to measure milk from all four quarters - if three quarters should measure 290 units and one 250 units, then there is a risk of infection in the lower reading quarter.
5. **Milk-fat content has an effect** on the magnitude of **readings**. Fatty milk gives above average readings.
6. **Illnesses other than mastitis, pharmaceutical procedures as well as heat may influence readings**. The recommended interpretations given in the instruction manual are not applicable in such cases.
7. **Dirty, fat-coated electrodes negatively influence readings**. The first step to take when higher than normal results occur, is to thoroughly clean the electrodes of any fat.

8. **Cows with visible clinical mastitis** will generally return to a reading above 300 units because of the essential changes occurring in severely infected milk.
9. Although **frequent measurements increase reliability of the test**, it is sufficient for large herds to be examined at established intervals. The data obtained from these measurements provide the basis for selecting the infection-threatened cows, which should be thoroughly investigated.
10. **In herds where a large number of cows have milk samples below 300** units a mechanical check of the milking machine by a machinery specialist is desirable, as wide-spread lower readings throughout the herd can indicate over-milking of cows or acidosis is occurring.
11. Although there is the **correlation between somatic cell count** and the Mastitis Detector's readings, the Detector cannot be used for somatic cells counting.  
The somatic cell count of milk is an insufficient diagnostic tool for determining mastitis. Udder health disorders may be present even with a low somatic cell count. There are microbes such as nonhaemolytic Staphylococci, which cause only minimal increases in the somatic cell count and as they multiply they cause inflammation of the udder.

## CARE & CLEANING

### **WHEN USED REPETITIVELY, RINSING BETWEEN QUARTERS OR BETWEEN COWS IS NOT USUALLY NECESSARY.**

Milk does not wet or stick to the polypropylene cup and any minute droplets remaining after a shake are too insignificant to materially affect the next reading.

Rinsing with water to prevent milk-fat drying on the terminals, is only necessary when the detector is not in constant use. It is also recommended whenever subclinical mastitis has been detected in a milk sample to avoid affecting the next reading.

When not in use clean the Detector with methylated spirits or detergent, rinse and dry, and store in its box in a dry place at normal room temperature.

In regions with high iron, calcium or phosphate levels in water, the terminals in the measurement cup should be periodically cleaned with a mildly abrasive kitchen detergent to remove any chemical residue or oxidation.

The DRAMIŃSKI Electronic Mastitis Detector is a highly sensitive electronic instrument designed to measure very small changes in milk electrical resistance accurately. The instrument has been designed for use in the often hostile environment of the milking area and is sufficiently robust, shockproof, and waterproof to resist most normal working conditions.

However, as a sensitive electronic instrument, its accuracy is best maintained by keeping it clean and free from injury, physical abuse or extreme temperatures.

**DO NOT USE BOILING OR VERY HOT WATER** for washing or immersing the instrument; **DO NOT DRY THE DETECTOR IN AN OVEN OR MICROWAVE, OR STORE IN AN ENVIRONMENT ABOVE 50 DEGREES C AS DAMAGE TO ELECTRONIC COMPONENTS CAN RESULT.**

## BATTERY REPLACEMENT

Under normal use, power consumption is very low and a 9 volt alkaline battery will provide many thousands of tests.

Lower battery power is indicated by the icon flashing on the LCD. A new battery is required when the icon starts flashing.

To change the battery:

1. Unscrew both screws which fix the small cover on the end of the handle.
2. Remove the battery from inside the handle and unclip the battery harness.
3. Clip the harness onto a new 9 volt or 9 volt alkaline battery, insert the battery into the handle.
4. Replace the cover, ensure the rubber gasket is in place and tighten the screws.

## TECHNICAL DATA

Total weight:	approx. 300 gm
Power:	9 volt battery, standard or alkaline
Power consumption:	approx. 14 mA
Display:	LCD 4 blocks of 3 digits
Smallest unit:	10 units
Measurement range:	10 to 990 units
Storage and work temperature:	0 to 50° C
Max. humidity:	95%
Low - battery indicator:	blinking battery icon

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